

22

7,9

42,0

9,0

▲	L	I
		☾

12

7,8

43,0

10,0

▲	L	I
		☾

18A

8,0

44,0

9,6

▲	L	I
		☾

5

7,7

42,5

10,0

▲	L	I
		☾

8A

8,0

44,0

10,0

▲	L	☾
		☾

17

7,5

44,5

10,5

▲	L	I
		☾

16

8,0

45,0

10,0

▲	L	I
		☾


19A

9,0

47,0

10,7

▲	L	☾
		☾



2

7,2

42,0

9,5

		☾
●		☾
P		

10

8,0

42,0

9,5

		I
●		☾
P		

33A

7,8

43,0

8,5

		☾
●		☾
P		

24

7,8

45,0

9,0

		☾
●		☾
P		

38A

8,0

44,0

10,0

		☾
●		☾
P		

36

8,3

45,0

9,0

		☾
●		☾
P		

37A

8,5

45,5

9,8

		I
●		☾
P		

30A






8,3



45,0

10,5




		☾
●		☾
P		








FLAT-PLANE



*Major Dent*




SHADE GUIDE  
2C-3R  
A1-D4




**MA**

**MEDIUM-ANATOMIC**




**ARTICULATION COMBINATION CHART**  
**ARTIKULATIONS-TABELLE**  
**TABEAU D'ARTICULATION**  
**TABLA DE ARTICULACION**  
**TABELLA DI ARTICOLAZIONE**  
**OVERSIGHT OVER ARTIKULATION**  
**PURENTATAULUKKO**  
**ΠΙΝΑΚΑΣ ΑΡΘΡΩΣΕΩΝ**  
**ARTICULATIETABEL**  
**TABELA DE ARTICULAÇÃO**  
**TABELL ÖVER LEDRÖRELSE**


			<b>FP</b>	<b>MA</b>
			FP (Flat-Plane)	MA (Medium-Anatomic)
2	22	3A	75R	52 59
3A	22	3A	74R	52 —
5	5	—	74R	56 65
7A	7A	—	74R	55A 60A
8A	8A	36	75R	55A 59
10	22	10	54A	56 60A
12	7A	—	75R	— 65
16	27A	5	75R	68A 59
17	5	36	75R	55A 60A
18A	7A	—	74R	68A 59
19A	19A	16	75R	68A —
22	22	—	54A	56 60A
24	36	37A	75R	68A 65
27A	27A	37A	74R	68A 59
29A	27A	36	75R	68A 59
30A	16	40A	75R	68A 59
33A	10	7A	74R	52 59
35A	36	37A	74R	55A 60A
36	27A	36	75R	68A 65
37A	37A	—	74R	68A 65
38A	5	8A	75R	68A 59
39	5	36	74R	— 59
40A	19A	40A	75R	68A —
45A	19A	—	75R	68A —




mm	mm			
46/49	6,90/7,60	5-17	2-33A	3A-35A
50/52	7,60/8,00	5-17	10-24	35A
53/54	8,00/8,30	18A-8A	10-38A-36	—
55/56	8,30/8,60	16	36-30A	27A
57/61	8,60/9,40	19A	37A	39-29A-40A



mm	mm			
30/31	37,5/40,7	—	—	—
32	40,4/43,2	22-12-5	2-10-3A	3A
33	42,2/44,5	12-18A-5-8A-17	33A-38A	35A-7A
34/35	44,0/47,2	18A-8A-17-16-19A	24-38A-36-37A-30A	35A-7A-27A-39-29
36/37	46,7/49,8	19A	—	27A-39-40A-45A



mm	mm			
4,5/4,7	6,70/7,60	5-17	2-33A	3A-35A
4,8/4,9	7,20/8,00	22-12-18A-5-17	2-10-33A-24	3A-35A-7A
5,0/5,1	7,50/8,30	22-12-18A-5-8A-17	10-33A-24-38A	3A-35A-7A
5,2/5,3	7,80/8,50	12-18A-8A-16	10-24-38A-36-30A	7A-27A
5,4/5,7	8,30/9,30	16-18A	36-37A-35A	27A-39-29A-40A



mm	mm			
111/120	36,7/40,0	—	—	—
120/125	40,0/42,0	22	2	3A
126/130	42,0/43,3	22-5	22-33A-10	—
130/135	43,3/45,3	18A-8A-17-16	33A-24-38A-36	35A-7A
136/140	45,3/46,7	—	37A	27A-29A
140/145	46,7/48,7	19A	—	39-40A
146/155	48,7/52,0	—	—	45A

CHOICE SYSTEM IN RELATION TO THE SIZE ■ AUSWAHLSYSTEM NACH DEM KÖRPERBAU ■ SYSTEME DE CHOIX EN RELATION A LA TAILLE ■ METODO DE ELECCION EN BASE A LA COMPLEXION FISICA ■ MISURAZIONE IN BASE ALLA CORPORATURA ■ MÅLING EFTER KROPPSBYGGNINGEN ■ MITOITUS RUUMIINRAKENTEEN ■ PERUSTEELLA ■ ΜΕΤΡΗΣΗ ΒΑΣΕΙ ΣΩΜΑΤΙΚΗΣ ΔΙΑΠΛΑΣΗΣ ■ METINGEN OP BASIS VAN LICHAAMSSPECIFIEKE VORMEN ■ ΜΕΔΙCΑΟ ΣΕΓΥΝΔΟ Α ΕΣΤΑΤΥΡΑ ■ MÄTNING PÅ BASIS AV KROPPSBYGGNAD



	Major Dent Flat-Plane	Major Dent Medium-Anatomic	Major Plus High-Anatomic	Major Plus Technical
Koko- ja osaproteesi - Ολική και μερική τεχνητή οδοντοστοιχία - Totale en gedeeltelijke prothese - Prótese total e parcial - Hel- och delprotiser	<b>FP</b> ●	<b>MA</b> ●	<b>HA</b> ●	<b>T</b> ●
Yhden leuan kokoproteesi - Ολική τεχνητή οδοντοστοιχία με φυσικό ανταγωνιστή - Totale prothese met natuurlijke antagonist - Prótese total com antagonista natural - Helprotiser med naturlig antagonist		<b>MA</b> ●	<b>HA</b> ●	
Ranka- ja yhdistelmäproteesi - Τεχνητή οδοντοστοιχία σκελετική και συνδυασμένη - Uitneembare of gecombineerde prothese - Prótese esquelética e combinada - Skelettproteser och kombinerade proteser		<b>MA</b> ●	<b>HA</b> ●	
Τεχνητή οδοντοστοιχία προσωρινή - Tijdelijke prothese - Prótese provisória - Provisoriska proteser	<b>FP</b> ●	<b>MA</b> ●	<b>HA</b> ●	<b>T</b> ●
Τεχνητή οδοντοστοιχία επί εγκατάστασης μεταλλικής λάμας - Prothese met implantaat - Prótese sobre estrutura - Proteser på fästen	<b>FP</b> ●	<b>MA</b> ●	<b>HA</b> ●	<b>T</b> ●
Proteesi, jossa lingualisesti asetellut hampaat - Τεχνητή οδοντοστοιχία με μοντάρισμα lingualized - Prothese met linguale bevestiging - Prótese com montagem lingualizado - Proteser med lingualmontering		<b>MA</b> ●	<b>HA</b> ●	