

MATERIAL SAFETY DATA SHEET – Edition 01  
28<sup>th</sup> April 2009

1. Identification of the product and of the company:-

Identification of the product:  
Non woven swabs sterile and non sterile

Company/undertaking identification:  
SYLAMED  
3 square Maubeuge  
75009 Paris - France

2. Composition of the product:-

70% viscose – 30% polyester.

3. Hazards Identification of the product:-

None if used for the intended purpose.

4. First Aid Measures:-

After Inhalation:  
Not relevant

After Skin Contact:  
Not relevant

After Eye Contact:  
Not relevant

After Ingestion:  
Not relevant

5. Fire-fighting Measures:-

Material is flammable and suitable measures should be taken to avoid any naked flames near the material.

Suitable extinguisher should be used for fighting fire.

6. Accidental release measures:-

Personal precautions:  
Not required

Environmental precautions:  
Not required

Methods of cleaning up/of removing  
Should be disposed of sensibly along with all other combustible materials like paper and plastic.

7. Handling and Storage:-

Handling:  
No special handling required

Storage:  
Should be kept clean, in a dry area and away from excessive moisture.

8. Exposure Controls / Personal Protection

No special actions required.

9. Physical and Chemical Properties:-

**Standards : Pharmacopée européenne II - Norme européenne EN 1644-1 et EN 1644-2**

**-Physiques**

Poids	40 g / m <sup>2</sup> ± 2,5 g
Dimensions pliée	± 2,5 %
Nombre d'épaisseurs	4
Blancheur	> 80 A°
Capacité d'absorption	10 g/g
Temps d'absorption	≤ 10 s
Élongation	MD = 115 ± 10 % CD = 35 ± 20 %

**-Chimiques**

PH	4,4 ≤ pH ≤ 8,2
Fluorescence	absence
Substances tensioactives	≤ 2 mm
Amidon / dextrine	absence
Substances hydrosolubles	≤ 0,5 %
Substances solubles dans l'éther	≤ 0,5 %
Cendres sulfuriques	≤ 2 %

10. Stability and Reactivity:-

None if used for intended purpose.

11. Toxicological Information:-

None Known

12. Ecological Information:-

Product is not bio-degradable.

13. Disposal Considerations:-

Waste incineration with the approval of the local authority responsible.

14. Regulatory Information:-

Not applicable.

15. Other Information:-

This information is based on our current level of knowledge and relates to the product in the state in which it is delivered. It is intended to describe our products from both the point of view of safety requirements and is not intended to guarantee any particular properties.